

LAB: SNC 2DI- Properties of Strong and Weak Acids and Bases

	A		B		C		D	
Litmus Paper • add ~1 cm strip • record colour and whether it is an acid/base	Red	Blue	Red	Blue	Red	Blue	Red	Blue
pH Paper • add ~1 cm strip • record colour and corresponding pH								
Phenolphthalein Indicator • add 2 drops to 5 mL of solution in small test tube • record colour and whether it is an acid/base								
Bromthymol Blue • add 2 drops to 5 mL of solution in small test tube • record colour and whether it is an acid/base								
Universal Indicator • add 2 drops to 5 mL of solution in small test tube • record colour and corresponding pH								

QUESTIONS:

- What does pH actually measure?
- What is the advantage of using Universal indicator or pH paper over litmus papers and phenolphthalein?
- The following mystery liquids were tested with various indicators. Classify them as acids, bases or neutral:
 - Red litmus paper turns blue, phenolphthalein turns pink
 - Bromthymol blue turns blue, blue litmus paper turns blue
 - Blue litmus paper turns blue, phenolphthalein stays clear
 - Phenolphthalein stays clear, bromthymol blue turns yellow
 - Red litmus paper turns red, phenolphthalein stays clear
 - Blue litmus paper turns red, phenolphthalein stays clear