

DETERMINATION OF RATE LAW EXPRESSIONS

Purpose: To determine the rate law with respect to concentration of hydrochloric acid when reacted with magnesium.

Materials: stop watch or clock 6 molar HCl
flasks graduated cylinder
4 three centimetre pieces of magnesium ribbon

Pre-lab: Do the calculations to produce the proper HCl solutions.

Procedure: 1. Prepare 50mL of 0.5M, 1M, 1.5M and 2M HCl

2. Place 50mL of 2M HCl into a 250 mL erlenmeyer flask.
(volume of flask is not critical)

3. Noting the time, drop a 3 cm piece of Mg into the 2M HCl

4. Swirl and/or shake the flask to obtain uniform results.

5. Record the time when all the ribbon has reacted.

6. Repeat steps 2 to 5 for each concentration of HCl

7. Calculate the rate of each reaction as $1/t$. There is no specific property being measured but time is the denominator.

8. Plot the rate of reaction ($1/t$) vs concentration of HCl. ($1/t$ on y-axis, conc. on x-axis)

9. Plot the rate of reaction ($1/t$) vs concentration² of HCl. ($1/t$ on y-axis, conc². on x-axis)

10. Plot the rate of reaction ($1/t$) vs concentration³ of HCl. ($1/t$ on y-axis, conc³. on x-axis)

11. The straightest line is the direct correlation between rate and conc^x.
Using this line state the rate law for the reaction. ie. Rate = $k [\text{HCl}]^x$

12. Calculate k for the reaction including units.