

## 8.4 Acid-Base Titration

- In an acid-base titration, the concentration of an acid or base solution of unknown concentration is determined by the delivery (from a buret) of a measured volume of a solution of known concentration (the titrant). If the sample in the flask is an acid, the titrant used is a base, and vice versa.
- Titration: the precise addition of a solution in a buret into a measured volume of a sample solution.
- Titrant: the solution in a buret during the titration.
- Sample: the solution being analyzed in a titration.
- Accuracy is essential. The titrant needs to be standardized. Primary standards are used to measure a titrant's concentration. Sodium carbonate ( $\text{Na}_2\text{CO}_{3(s)}$ ) is a base used to standardize acids. Potassium hydrogen phthalate ( $\text{KHC}_7\text{H}_4\text{O}_{4(s)}$ ) is an acid used to standardize bases.
- Primary Standard: a chemical available in a pure and stable form, for which an accurate concentration can be prepared; the solution is then used to determine precisely, by means of titrating, the concentration of the titrant.
- Equivalence point: in a titration, the measured quantity of titrant recorded at the point at which chemically equivalent amounts have reacted.
- Endpoint: the point in a titration at which a sharp change in a measurable and characteristic property occurs; e.g. a colour change in an acid-base indicator.
- How to perform a basic titration where the titrant is a dilute strong base and the sample is a dilute strong base.
  - Create your stock titrant solution following good lab practices to minimize errors.
  - Rinse the buret with the stock solution and allow it to flow through the wide-open nozzle to knock out any air bubbles.
  - Fill the buret, keeping in mind that you can waste a great deal of time trying to get it to zero, therefore don't. Just fill to near the top and record your starting volume.
  - Standardize your titrant by accurately measuring and creating a primary standard solution of known concentration.
  - Add an indicator that will change colour at the equivalence point. The endpoint of the indicator should be close the equivalence point.

- Carefully add the titrant to the primary standard until you reach the endpoint. Record the volume. Repeat until results are consistent ( $\pm 0.2$  mL).
- Calculate the amount of titrant delivered.
- Calculate the concentration of the titrant.
- Perform the titration on the unknown and repeat until results are consistent.
- Calculate the concentration of the unknown.

### **Titrating a Strong Acid with a Strong Base**

- Consider the neutralization of HCl with NaOH:
  - $\text{HCl}_{(\text{aq})} + \text{NaOH}_{(\text{aq})} \rightarrow \text{H}_2\text{O}_{(\text{l})} + \text{NaCl}_{(\text{aq})}$       molecular equation
  - $\text{H}_3\text{O}^+ + \text{OH}^- \rightarrow 2\text{H}_2\text{O}_{(\text{l})}$       net ionic equation
  - $\text{H}^+_{(\text{aq})} + \text{OH}^- \rightarrow \text{H}_2\text{O}_{(\text{l})}$       abbreviated net ionic eq.
- E.g. In a titration, 20.00 mL of 0.300 mol/L  $\text{HCl}_{(\text{aq})}$  is titrated with standardized 0.300 mol/L  $\text{NaOH}_{(\text{aq})}$ . What is the amount of unreacted  $\text{HCl}_{(\text{aq})}$  and the pH of the solution after the following volumes of  $\text{NaOH}_{(\text{aq})}$  have been added?
  - (a) 0 mL
  - (b) 10.0 mL
  - (c) 20.0 mL

*Strong acid and strong base, therefore, the ions of interest are  $\text{H}^+$ ,  $\text{OH}^-$ , and  $\text{H}_2\text{O}$ . We can use the abbreviated net ionic equation to represent the reaction.  $\text{H}^+_{(\text{aq})} + \text{OH}^- \rightarrow \text{H}_2\text{O}_{(\text{l})}$*

- (a) No NaOH has been added therefore we are only concerned with HCl. Since it is a strong acid it will completely ionize and we know we will have 0.300 mol/L of  $\text{H}^+$  and 0.300 mol/L of  $\text{Cl}^-$

$$[\text{H}^+] = 0.300 \text{ mol/L}$$

$$\text{pH} = -\log [\text{H}^+] = -\log [0.300] = 0.5$$

$$n = v \times C = 0.020 \text{ ml} \times 0.300 \text{ mol/L} = 0.006 \text{ mol}$$

*Before we add any titrant the pH is 0.5 and we have 0.006 moles of HCl (but since it completely disassociates we have 0.006 mol of  $H^+$  and 0.006 mol of  $Cl^-$ ).*

- (b) By adding 10 mL of 0.300 mol/L  $NaOH_{(aq)}$  to 20.00 mL of 0.300 mol/L  $HCl_{(aq)}$  we can deduce that the reaction has reached the halfway point. We need to figure out how much of the acid was neutralized.

NaOH

$$n = v \times C = 0.010 \text{ L} \times 0.300 \text{ mol/L} = 0.003 \text{ mol}$$

Therefore we have 0.003 mol of NaOH but since it completely disassociates we have 0.003 mol of  $Na^+$  and 0.003 mol of  $OH^-$ .

Using the equation,  $H^+_{(aq)} + OH^- \rightarrow H_2O_{(l)}$ , we know we started with 0.006 mol of  $H^+$  and we added 0.003 mol of  $OH^-$ . Therefore the  $H^+$  will combine with  $OH^-$  using up all the  $OH^-$  and leaving behind 0.003 mol of  $H^+$ .

We must also take into consideration the new volume: 20 mL + 10 mL = 30 mL

$$[H^+] = 0.003 \text{ mol} / 0.03 \text{ L} = 0.100 \text{ mol/L}$$

$$pH = -\log [H^+] = -\log [0.100] = 1.0$$

*We have 0.003 mol of  $H^+$  after adding 10 mL of NaOH and it will have a pH of 1.0.*

- (c) By adding 20 mL of 0.300 mol/L  $NaOH_{(aq)}$  to 20 mL of 0.300 mol/L  $HCl_{(aq)}$  we can deduce that the reaction has reached the endpoint. All the  $OH^-$  and  $H^+$  will have combined to form water and since  $Na^+$  and  $Cl^-$  do not have an influence on pH we can determine that  $[H^+] = 1.0 \times 10^{-7}$  (from the auto-ionization of water) and  $pH = 7$ .

- Look at Figure 2 on page 599 to see the titration curve for this question.
- The midpoint will have a rapid change in pH.

## Titrating a Weak Acid with a Strong Base

- These are the most challenging type of questions. Since we are also working with a weak acid we have to determine how equilibrium will exert itself after the base is added.
- E.g. In a titration of 20.00 mL of 0.300 mol/L  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  (acetic acid) with standardized 0.300 mol/L  $\text{NaOH}_{(\text{aq})}$ , what is the amount of unreacted  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  and the pH of the solution:
  - (a) before titration begins
  - (b) after 10.00 mL of NaOH is added
  - (c) after 20 mL of NaOH is added.

*We need to remember that this is a weak acid and will not disassociate easily and equilibrium will re-exert itself.*

(a) Before the titration only  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  has to be considered which makes it a simple equilibrium problem.

	$\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$	$\leftrightarrow$	$\text{H}^+_{(\text{aq})} +$	$\text{C}_2\text{H}_3\text{O}_2^-_{(\text{aq})}$
Initial	0.300		0	0
Change	-x		+x	+x
Equilibrium	0.300-x		x	x

$$K_a = 1.8 \times 10^{-5} = \frac{[\text{H}^+_{(\text{aq})}][\text{C}_2\text{H}_3\text{O}_2^-_{(\text{aq})}]}{[\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}]} = \frac{(x)(x)}{0.300 - x}$$

*We can use the hundred rule and assume that  $0.300 - x = 0.300$*

$$\frac{x^2}{0.300} \cong 1.8 \times 10^{-5}$$

$$x^2 \cong 5.4 \times 10^{-6}$$

$$x \cong 2.3 \times 10^{-3}$$

$$[\text{H}^+] = 2.3 \times 10^{-3} \text{ mol/L}$$

$$\text{pH} = -\log [\text{H}^+] = -\log [2.3 \times 10^{-3}] = 2.62$$


$$n = v \times C = 0.020 \text{ ml} \times 0.300 \text{ mol/L} = 0.006 \text{ mol}$$

*The starting pH is 2.62 and the starting amount is 0.006 mol.*

(b) After adding 10 mL of NaOH.

*Perform the stoichiometric calculations first.*

The  $\text{OH}^-_{(\text{aq})}$  will react with  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  to form acetate and water.



How much  $\text{OH}^-$  reacts with the acid (assume that NaOH disassociates completely)?

$$n = V \times C = 0.01 \text{ L} \times 0.300 \text{ mol/L} = 0.003 \text{ mol}$$

We have already calculated that we have 0.006 moles of  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$ , therefore 0.003 moles of  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  will be neutralized by the  $\text{OH}^-_{(\text{aq})}$ .

$0.006 \text{ mol} - 0.003 \text{ mol} = 0.003 \text{ mol}$  (*All the  $\text{OH}^-_{(\text{aq})}$  is consumed and 0.300 mol of  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  remain.*)

$20 \text{ mL} + 10 \text{ mL} = 30 \text{ mL}$  (*the new volume*)

*Perform equilibrium calculations*

Before equilibrium exerts itself we know a few things. We have 0.003 mol of  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  and 0.003 mol of  $\text{C}_2\text{H}_3\text{O}_2^-_{(\text{aq})}$  in 0.03 L to give 0.100 mol/L of each. We now need to figure out how equilibrium will exert itself.

	$\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$	$\leftrightarrow$	$\text{H}^+_{(\text{aq})} +$	$\text{C}_2\text{H}_3\text{O}_2^-_{(\text{aq})}$
Initial	0.100		0	0.100
Change	-x		+x	+x
Equilibrium	$0.100 - x$		x	$0.100 + x$

$$K_a = 1.8 \times 10^{-5} = \frac{[H^+_{(\text{aq})}][C_2H_3O_2^-_{(\text{aq})}]}{[HC_2H_3O_{2(\text{aq})}]} = \frac{(x)(0.100 + x)}{0.100 - x}$$

*We can use the hundred rule and assume that  $0.100 - x = 0.100$  and  $0.100 + x = 0.100$*

$$\frac{(x)(0.100)}{0.100} \cong 1.8 \times 10^{-5}$$

$$x \cong 1.8 \times 10^{-5}$$

$$[\text{H}^+] = 1.8 \times 10^{-5} \text{ mol/L}$$


$$\text{pH} = -\log [\text{H}^+] = -\log [1.8 \times 10^{-5}] = 4.74$$

*After 10 mL of NaOH, the pH is 4.74 and the amount of acid is 0.003 mol.*

(c) After adding 20 mL of NaOH at the equivalence point.

*Perform the stoichiometric calculations first.*

The  $\text{OH}^-_{(\text{aq})}$  will react with  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  to form acetate and water.



How much  $\text{OH}^-$  reacts with the acid (assume that NaOH disassociates completely)?

$$n = V \times C = 0.02 \text{ L} \times 0.300 \text{ mol/L} = 0.006 \text{ mol}$$

We have already calculated that we have 0.006 moles of  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$ , therefore all of  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  will be neutralized by the  $\text{OH}^-_{(\text{aq})}$ .

$0.006 \text{ mol} - 0.006 \text{ mol} = 0 \text{ mol}$  (*All the  $\text{OH}^-_{(\text{aq})}$  is consumed and none of the  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  remain.*)

$20 \text{ mL} + 20 \text{ mL} = 40 \text{ mL}$  (*the new volume*) However, we will have all the  $\text{HC}_2\text{H}_3\text{O}_{2(\text{aq})}$  converted into  $\text{C}_2\text{H}_3\text{O}_2^-_{(\text{aq})}$ .

*Perform equilibrium calculations*

Before equilibrium exerts itself we know a few things. We have 0.006 mol of  $\text{C}_2\text{H}_3\text{O}_2^-_{(\text{aq})}$  in 0.04 L to give 0.150 mol/L. We now need to

figure out how equilibrium will exert itself. But the reaction now changes. We have 100% acetate and we have to figure out the equilibrium for it.

	$C_2H_3O_2^-$ (aq)	+	$H_2O$ (l)	$\leftrightarrow$	$HC_2H_3O_2$ (aq)	+	$OH^-$ (aq)
Initial	0.150		*		0		0
Change	-x		*		+x		+x
Equilibrium	0.150-x		*		x		x

$$K_b = \frac{K_w}{K_a} = \frac{1.0 \times 10^{-14}}{1.8 \times 10^{-5}} = 5.6 \times 10^{-10}$$

$$K_b = 5.6 \times 10^{-10} = \frac{[HC_2H_3O_2] [OH^-]}{[C_2H_3O_2^-]} = \frac{(x)(x)}{0.150 - x}$$

*We can use the hundred rule and assume that  $0.150 - x = 0.150$*

$$\frac{x^2}{0.150} \cong 5.6 \times 10^{-10}$$

$$x \cong 9.2 \times 10^{-6}$$

$$[OH^-] = 9.2 \times 10^{-6} \text{ mol/L}$$

$$pOH = -\log [OH^-] = -\log [9.2 \times 10^{-6}] = 5.03$$

$$pH = 14 - 5.03 = 8.97$$

*After 20 mL of NaOH, the pH is 8.97 and the amount of acid is 0 mol.*

### **Titrating a Weak Base with a Strong Acid**

- The questions are similar to the one above.

### **Acid-Base Indicators**

- Indicators will change colour at a specific pH.
- Transition Point: the pH at which an indicator changes colour

$$K_{In} = \frac{[In^-] [H_3O^+]}{HIn}$$

- See table 7 on page 609 for information about the indicators.
- You can calculate which indicator will be the best to use by calculating the pH at a reactions equivalence point.

### **Polyprotic Acid Titrations**

- These titrations are the same as monoprotic titrations except there are multiple endpoints.

### **Homework**

- Practice 1,2,3,4,5,6,10
- Questions 1,2,3,4,5,6,7,8,9,10